



## Changes to the Australian Hearing FM range, June 2009.

Australian Hearing has made several changes to its range of Top Up FM systems (those that require a client contribution) due to changes in Phonak's FM range.

Two products have been removed from the paediatric FM Top Up range

1. Phonak iLink combined hearing aid & FM receiver.
2. EduLink ear level receiver for unaided clients has been discontinued and will be replaced by the iSense Micro.

Clients who have previously purchased EduLink or iLink:

The EduLink or iLink devices that are within warranty will be covered by the warranty. Service parts to maintain the devices will be available for **5 years** after the phase out date for both EduLink and iLink.

**Additions:** Two new top up products have been added to the range that is available to children who do not wear hearing aids but who have a mild, unilateral or fluctuating conductive hearing loss.

Safety: Both these FM receivers contain small parts and are not suitable for children under 3 years of age.

The iSense products shown below offer Dynamic FM features when used in conjunction with the Inspiro transmitter.

### **iSense Micro.**

iSense Micro replaces EduLink as the ear level receiver option for unaided children who will benefit from an FM system. The receivers are available in 6 colours (see below)



iSense Micro does not have an on-unit volume control. However a volume control can be activated on the Inspiro transmitter.

### **3. iSense Classic**

iSense Classic is a lightweight option for children who wish to use an FM receiver with earphones. As shown in the picture below, the Classic uses the same housing as the MyLink receiver; *however the two are not the same device*. The iSense Classic has an inbuilt, rechargeable battery. There are no colour options.

